



Recipient Name: _____


Amtryke Assessment Form

(Must be filled out completely by therapist)

Age: _____ Weight (lbs.):* _____ Height (inches):* _____
Diagnosis(es):* _____

*This private information is only used to help appropriately fit the rider.

RIDER'S MEASUREMENTS				
Arm Measurements (inches)				Total Length
Left	A to B:		B to C:	
Right	A to B:		B to C:	
Trunk	A to D:			
Leg Measurements (inches)				Total Length
Left	D to E:		E to F:	
Right	D to E:		E to F:	



A	Center of Shoulder
B	Center of Elbow
C	Center of Digit Crease
D	Center of Hip
E	Center of Knee
F	Bottom of Foot

Arm Length & Leg Length Measurements are critical to correct Amtryke Selection Sizing Chart is available online:
www.ambucs.org/riders/wish-list/sizing-chart/

Notes on Provided Measurements (if any):

Helmet Sizing	
Size	Measurement (head circumference)
XXS	18.5" to 19.5"
XS	20.5" to 22"
S/M	22" to 23.6"
L/XL	23.6" to 35.75"

Therapist Name: _____ Is this the treating therapist? Yes No

Credentials: _____

Phone: _____ Email: _____

Facility Name: _____

Street Address: _____ City: _____ State: _____ Zip: _____

Is this facility an Amtryke Evaluation Site? Yes No Not Sure

Therapist comments concerning recipient or goals:

By signing below, you are signifying that in your professional opinion this rider would benefit from an Amtryke. You assume no liability.

Therapist Signature: _____ Date: _____

Ship Amtryke To

Name/Facility: _____	Phone: _____
Street Address: _____	
City: _____	State: _____ Zip: _____



Recipient Name: _____

Thanks for choosing an Amtryke therapeutic tricycle!

In order to accommodate the widest variety of people and disabilities, Amtryke offers many tryke models and each can be customized in a variety of ways. Following the steps below will help you choose to the perfect tryke for your client from what might seem like a dizzying array of options.

Step 1: Fill out the first page of the Amtryke Assessment Form.

Step 2: Choose the way the tryke will be propelled: **Hand & Foot, Foot, or Hand.** Your choice should be based on the rider's ability and therapy goals.

Hand & Foot trykes improve coordination, strength and range of motion. Using all four extremities helps with weakness in any area, even general weakness, and can positively affect tone.

Foot trykes were developed in response to requests from therapists for a traditional tricycle for riders with special needs.

Hand trykes are designed for persons whose lower limbs lack function or those who need special therapy for the upper extremity.

Step 3: Take rider's measurements from the front of this Assessment Form and compare them to our Sizing Chart. (This will narrow the choices considerably)

Step 4: Choose a drive. (If it doesn't mention a choice, then ignore this step.)

Amtrykes come with two drive possibilities: **fixed** drive or **geared** drive. Tryke models have been carefully designed so the drive matches the therapeutic goals of the equipment; therefore all drives are not available on all trykes.

A **fixed** drive, commonly known as a 'fixie,' works on a mechanical level to help individuals make a full pedal rotation. The foot crank is constantly in motion for full therapeutic effect. Coasting is not possible with a positive drive; when limb motion stops, the bike does as well.

Riders of **geared** drive trykes must be able to make a complete pedal rotation on their own. This tryke is suitable for riders who need help with balance and a stable sitting position and who have the cognitive ability to steer successfully and apply the coaster brake or hand brake. The key feature of geared trykes is the ability to coast.

Step 5: Chose any adaptations and/or accessories needed by the rider. *Each tryke model can only be customized in the ways noted in its own model section in the Tryke Selection Forms or with the generic accessories listed below. If a customization option is not listed, it is because of design or other practical constraints in offering it on a particular model.*

GENERIC ACCESSORIES (not model specific)				
Safety Items	<input type="checkbox"/> Safety Flag	<input type="checkbox"/> NEW! Medical Cart (carries O ₂ bottle/ventilator/tubing)	<input type="checkbox"/> Water Bottle w Cage	<input type="checkbox"/> License Plate
Leg & Foot Items	<input type="checkbox"/> Foot Cups (pair): <input type="checkbox"/> Small <input type="checkbox"/> Medium	<input type="checkbox"/> Knee Adductor Strap: <input type="checkbox"/> Small <input type="checkbox"/> Large <input type="checkbox"/> Medium	<input type="checkbox"/> Pedal Block (1 = 3/4") _____ qty	
Hand Items	<input type="checkbox"/> Variable Range of Motion Kit (only for Hand & Foot Cycles)	<input type="checkbox"/> Wrist Wraps (Includes right & left) <input type="checkbox"/> X-Small <input type="checkbox"/> Large <input type="checkbox"/> Small <input type="checkbox"/> Medium	<input type="checkbox"/> Wrist Brace Mitt: <input type="checkbox"/> Right <input type="checkbox"/> Left <input type="checkbox"/> XX-Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Small <input type="checkbox"/> Small	



GRAPEVINE
1-877-AMTRYKE
grapevineambucs.org

AMBUCS

Recipient Name: _____

HAND-FOOT TRYKES

<input type="checkbox"/> New! AM-10 HAND-FOOT TRYKE (12" wheels & 3" crank arms) 50-HFC-0105 Arm length 12-16", leg length 14-21"			
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat			
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System			
Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest			
<input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**			
AM-10 Accessories:	<input type="checkbox"/> Rear Steering Kit	<input type="checkbox"/> Pull/Steering Bar	<input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
	<input type="checkbox"/> Push Bar	<input type="checkbox"/> Heavy-Duty Push Bar	<input type="checkbox"/> Cube Separator

<input type="checkbox"/> AM-12Small HAND-FOOT TRYKE (12" wheels & 3" crank arms) 50-HFC-0110 Arm length 12-16", leg length 16-21"			
Tryke comes standard with rear steering!			
<input type="checkbox"/> Option #1 – Standard Seating System: Bench Seat & 1600 Simple Seat Back			
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System			
Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest			
<input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**			
<input type="checkbox"/> Option #3 – Upgraded Seating System: Blue Bucket Seat			
AM-12S Accessories:	<input type="checkbox"/> Pull/Steering Bar	<input type="checkbox"/> Toe Pedal Leveler Pulley	<input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
		<input type="checkbox"/> Cube Separator	<input type="checkbox"/> Exercise pedals

<input type="checkbox"/> AM-12 HAND-FOOT TRYKE (12" wheels & 4" crank arms) 50-HFC-0210 Arm length 14-23", leg length 21-29"			
Tryke comes standard with rear steering!			
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back			
Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")			
<input type="checkbox"/> Option #2 – Upgraded Seating System: ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals)			
Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")			
Trunk Support Upgrades: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad			
<input type="checkbox"/> Option #3 – Upgraded Seating System: Gray Bucket Seat			
<input type="checkbox"/> Option #4 – Upgraded Seating System: Snappy Seat System			
Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest			
<input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back** <input type="checkbox"/> Cube Separator			
AM-12 Accessories:	<input type="checkbox"/> Pull/Steering Bar	<input type="checkbox"/> Toe Pedal Leveler Pulley	<input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
	<input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5"		<input type="checkbox"/> Exercise pedals

*cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

<input type="checkbox"/> AM-16 HAND-FOOT TRYKE (16" wheels & 5" crank arms) 50-HFC-0411 Arm length 18-27", leg length 24-36"			
Now comes standard with rear steering!			
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back			
Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")			
<input type="checkbox"/> Option #2 – Upgraded Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals)			
Seat Bottom Upgrades: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")			
Trunk Support Upgrades: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad			
<input type="checkbox"/> Option #3 – Upgraded Seating System: Gray Bucket Seat			
AM-16 Accessories:	<input type="checkbox"/> Pull/Steering Bar	<input type="checkbox"/> Toe Pedal Leveler Pulley	<input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
	<input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5"	<input type="checkbox"/> 9" Extender Tube	<input type="checkbox"/> Exercise pedals <input type="checkbox"/> XL Exercise pedals

Please mail, email or fax completed form to Grapevine Ambucs
Email: grapevineamtryke@gmail.com

Address: 2150 W. Northwest Highway Suite 114 #1089 Grapevine TX 76051
Fax: 855.827.1629

This Form, Amtryke Selection Form and the Request/Liability Waiver Form must be received by Grapevine Ambucs before placement is considered.



Recipient Name: _____

New! AM-20 HAND-FOOT TRYKE (Recumbent Single Speed) – 50-HFC-0610 Arm length 18-30", leg length 32-46"

Tryke is currently being used for research!

AM-20 Accessories:

Exercise Pedals

XL Exercise Pedals

Expanding Pedals

Wheelchair Seat Bar Ends

Swing Away Arms

Calf & Leg Supports (Expanding pedals)

Therapist Assembly Notes or Comments (if any):